

inLab Profire sinter furnace

Procedure for installing the software version R02H



Note:

The R02G update is suitable for all furnace versions of the inLab Profire series. Before installing the update, we recommend creating backup file of the furnace data on a USB stick. Please follow the instructions in the user manual /Section 5.6 "Data transfer via USB stick" and select the option "Backup data".

You have downloaded the software update R02H for the inLab Profire sinter furnace from the Internet site <https://www.dentsplysirona.com/inlab-downloads>.

The update provides the following:

- General bug fixes
- Improvements in the area of furnace control and user interface

The update is split in two steps:

1. Installation of the new software for the user interface
2. Installation of the new software for the control board

The installation will run automatically as soon as the update process was started.

Note:

Since a complex software package is installed, the update time is about 2.5 hours. The display freezes during the update process and does not show any update progress.

Please do not abort the update manually before the 2.5 hours have elapsed.

To install the software, please proceed as follows:

Copy the files **Build.zip** and **GS1_R02H.zip** contained in the folder "LAB-software-download-inLab-profired-R02H.zip" to the root directory of a USB stick.

Note:

Please make sure that the USB stick (at least USB 2.0) is formatted in FAT32, that there are no other files or folders on the USB stick and that the stick is not password protected.

Follow the instructions in the user manual/chapter 5.6 "Data transfer via USB stick" and select the option "Firmware update".

When the installation is complete, a message will appear indicating that the update has been successfully installed. Please restart the inLab Profire by using the power switch on the back of the furnace. Remove the USB stick and wait approx. 60 seconds before switching on the furnace again until the green LED on the front panel is no longer illuminated.

Please check the software version number in the furnace settings.

To do this, tap the "Settings" icon in the main menu and then tap on "Device data". Now "Software Version" "R02H/R02E.1" should be displayed.

If you have any questions, please contact your dealer or the Dentsply Sirona Customer Service.

Date: October 2020