

# inPost compatibility overview for scanning on multi-unit abutments from nt-trading

Implant manufacturer	Implant system	Platform	Multi-Unit Abutments	nt-trading: 2-CONnect	SIRONA	
			Series		inPost	REF.
<b>BioHorizons</b>						
	Tapered plus, Tapered internal, Tapered internal tissue level, Internal dental implant, Single stage dental implants	<span style="color: orange;">■</span> 3.5 <span style="color: green;">■</span> 4.5	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483
<b>Biomet 3i</b>						
	NanoTite Certain® PREVAIL, OSSEOTITE Certain PREVAIL, NanoTite Tapered Certain PREVAIL, NanoTite Parallel Walled Certain, OSSEOTITE Parallel Walled Certain, Full OSSEOTITE Parallel Walled Certain, OSSEOTITE XP Certain, Full OSSEOTITE XP Certain, NanoTite Tapered Certain, OSSEOTITE Tapered Certain, Full OSSEOTITE Tapered Certain	<span style="color: purple;">■</span> 3.4 <span style="color: blue;">■</span> 4.1	H-Serie	H 800 S H 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483
<b>Dentsply Implants (AstraTech, Friadent)</b>						
Astra Tech	OsseoSpeed TX	<span style="color: teal;">■</span> 3.5/4.0 <span style="color: purple;">■</span> 4.5/5.0	S-Serie	S 800 S S 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639
Friadent	Frialit / Xive	<span style="color: grey;">■</span> 3.4 <span style="color: orange;">■</span> 3.8 <span style="color: blue;">■</span> 4.5	T-Serie	T 800 S T 805 S T 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6559483 6551639 6551639
<b>MIS</b>						
	Seven (internal hex.)	<span style="color: purple;">■</span> Standard <span style="color: green;">■</span> Wide	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483
<b>Nobel Biocare</b>						
	<b>NobelReplace®</b> (Straight, Tapered, Platform Switching), <b>Replace Select</b> (Straight, Tapered, TC), <b>NobelSpeedy Replace</b>	<span style="color: orange;">■</span> RP <span style="color: blue;">■</span> WP	E-Serie	E 810 S E 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639
	Nobel Active, NobelReplace Conical Connection	<span style="color: red;">■</span> NP <span style="color: orange;">■</span> RP	F-Serie	F 800 S F 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639
<b>Osstem / Hiossen</b>						
	Osstem TS (USA: Hiossen ET)	<span style="color: green;">■</span> Standard USA:Regular	S-Serie	S 820 S	inPost for 2-CONnect KS61	6551639
<b>Straumann</b>						
	Standard (Roxolid), Standard Plus, Tapered Effect	RN (4,8 mm) WN (6,5 mm)	N-Serie	N 810 S N 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639
	Bone Level	NC (3,3 mm) RC (4,1/4,8 mm)	L-Serie	L 800 S L 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	6551639 6551639
<b>Zimmer</b>						
	Tapered Screw-Vent	<span style="color: green;">■</span> 3.5 <span style="color: purple;">■</span> 4.5	R-Serie	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	6559483 6559483

Technische Änderungen und Irrtümer vorbehalten, Best.-Nr. A91100-8945-01-7600, Printed in Germany, Dispo-Nr. 04605, 1015.Y0