

PROMARK®

ENDODONTIC MOTOR



THE MARK OF PERFORMANCE

Built upon a tradition of performance and reliability, the ProMark® is our most versatile and complete endodontic motor for rotary and reciprocation.

The brushless motor delivers the ample power and consistency of full A/C mains power and is compatible with a variety of e-type handpieces, including your Tulsa Dental Specialties 8:1 reduction Contra Angles (sold separately).

The ProMark is preloaded with the complete library of Tulsa Dental Specialties' rotary and reciprocation settings and is fully programmable for individual preferences. It can also be updated with future endodontic innovations for a longer return on your investment.

- ▶ Modern and attractive motor system
- ▶ Delivers smooth reciprocation
- ▶ Vivid 3-row LCD Display
- ▶ Integrated rotary & reciprocation file series library
- ▶ Auto-sensing global voltage compatibility (110/220 V)
- ▶ Well-weighted foot pedal with 360° activation
- ▶ Preset options allow five single file, file series, or custom file series
- ▶ Upgradable software & file library
- ▶ Handpiece calibration

Contra Angle sold separately.

DESCRIPTION	SKU
ProMark Endodontic Motor <i>Includes Console, Motor & Cord, 360° Foot Pedal, Handpiece Cradle & Bracket</i>	PMKEM1
ProMark Replacement Motor	PMKEM1MTR
ProMark Foot Pedal	PMKEM1PDL
Tulsa Dental Specialties 8:1 Mini-Head Push-Button Contra Angle	TUL-8M
Tulsa Dental Specialties 1:1 Contra Angle	TUL-1
1:5 Increaser Contra Angle (NSK)	AHP71

- ▶ Auto-Stop/Reverse/Forward in rotary mode
- ▶ Powerful 30,000 RPM brushless motor can be autoclaved
- ▶ With recommended 8:1 handpiece (sold separately):
 - ▶ Speed Settings: 50 to 3,600 RPM in rotary mode
 - ▶ Torque Settings: 30 to 1000 g-cm
- ▶ Capable of 150,000 RPM operation with 1:5* increaser handpiece
- ▶ Compatible with a broad range of e-type handpieces (1:1, 8:1, 16:1, 1:5*)
- ▶ 8:1 contra angle required for reciprocation

1-800-662-1202
For the latest information consult
www.TulsaDentalSpecialties.com

DENTSPLY
TULSA DENTAL
SPECIALTIES