

# Simplant®

## Simplant® Guide provides significantly improved accuracy

### Simplant® Guide


compared to freehand

### >3 times higher accuracy

**10** experienced\* surgeons placed implants freehand or with guides, in a model of the anterior maxilla. Accuracy of planned and actual implant placement was compared by scanning and analysis in Mimics Software.


Significantly higher accuracy evaluating:

**Angular Lateral Depth**



Actual and planned positions

Overall angular accuracy  
n= 40, mean, SD, min/max



**40**  
Guides\*\*

\* Experienced surgeons was defined as: minimum of 10 years of implant dentistry experience and 300 to 1000 implant installations yearly.

\*\* Tooth-supported Simplant guides with 4 mm high titanium inserts for either single or multiple edentulism.



Vermeulen J. The accuracy of implant placement by experienced surgeons: Guided vs freehand approach in a simulated plastic model. Int J Oral Maxillofac Implants 2016;E-pub Oct 14, doi: 10.11607/jomi.5065. [Abstract](#)



www.dentsplysirona.com

**Dentsply  
Sirona**  
Implants